

May 7, 2024

Test Update – Congenital CMV, Saliva (CCMVB)

Effective Date **Tuesday, May 21, 2024**

Description of Change The test methodology for the Congenital CMV Saliva Assay (test code: CCMVB) is undergoing enhancements to improve sensitivity and specificity. As part of this update, **the collection device will be upgraded from a dry swab to a swab submitted in Universal Transport Media (UTM).**

Collection Device

OLD DEVICE: FLOQswab in dry tube



NEW DEVICE: Universal Transport Media (UTM)



If you have additional questions regarding this test, please contact Technical Support between the hours of 8 a.m. and 4:30 p.m. at 877-402-4221.