

July 10, 2024

# TEST UPDATE

Vasoactive Intestinal Peptide		VIP			
Description of Change	Test code VIP (Vasoactive Intestinal Peptide) has been discontinued effective Monday, 7/1/24. Please order test code MISC for “Vasoactive Intestinal Polypeptide” or contact Technical Support at 877-402-4221 to discuss testing alternatives.				
Effective Date	<b>Immediately</b>				
Test Name	Vasoactive Intestinal Peptide				
Test Code	<table border="1"> <thead> <tr> <th>CURRENT TEST CODE</th> <th>INTERIM REPLACEMENT TEST CODE</th> </tr> </thead> <tbody> <tr> <td>VIP</td> <td>MISC</td> </tr> </tbody> </table>	CURRENT TEST CODE	INTERIM REPLACEMENT TEST CODE	VIP	MISC
CURRENT TEST CODE	INTERIM REPLACEMENT TEST CODE				
VIP	MISC				
CPT Code	84586				
Methodology	Direct Enzyme Immunoassay, Enzyme-Linked Immunosorbent Assay (ELISA)				
<b>Specimen Requirements</b>					
Minimum Volume	2 mL EDTA Plasma				
Container	Lavender Top EDTA Tube				
Collection	Collect specimen in lavender top EDTA tube. Centrifuge and transfer EDTA plasma to aliquot tube and immediately freeze.				
Special Instructions	Patient should be fasting 10-12 hours prior to collection of specimen. Patient should not be on any antacid medication or medications that affect intestinal motility for at least 48 hours prior to collection.				
Reference Range	Up to 36 pg/mL				

### Additional Information

If you have additional questions regarding these tests, please contact Technical Support between the hours of 8 a.m. and 4:30 p.m. For general questions, Customer Care is available to assist at any time.

### Technical Support/Customer Care

Call 877-402-4221

0062824\_v1